



Mizuho Head Frame Sterilization Case



- Designed for storage and sterilization of Sugita Multi-Purpose Head Frame or Radiolucent Head Frame
- Accommodates the entire retraction system
- Keeps system well protected and organized
- Extremely durable



Product Description



The two tier system accommodates the entire Radiolucent Head Frame and everything but the table attachment for the Sugita Multi-Purpose Head Frame. The base of the case is constructed from durable aluminum to protect the head frame during cleaning and sterilization. The upper tray contains an inner container to hold scalp hooks and separate compartments for the various accessories. The case for the Radiolucent Head Frame has a warning marked on it stating not to autoclave the device.

Ordering Information

Catalog Number

Sugita Multi-Purpose Head Frame Case

07-980-00

Radiolucent Head Frame Case

07-980-00R

Mizuho America, Inc.
133 Brimbal Ave.
Beverly, MA 01915
Tel: 800-699-CLIP (2547) or
978-921-1718
Fax: 978-921-4003
EMAIL: mizuho@mizuho.com
www.mizuho.com



© Mizuho is a registered trademark of Mizuho Ikakogyo Co., Ltd. Tokyo, 2002

©2002 Mizuho America, Inc. HFC001 Rev. 1 5K 06/02MW