

Micro Scissors for Cerebrovascular Applications

Precision-designed to allow cutting through thick vascular and scar tissue.
Serrated blades prevent tissue slippage and stay sharper longer.



Fukushima/Giannotta Micro Scissors

Traditional Micro Scissors

- Ideal for cutting thick vascular and scar tissue
- Round or flat handle
- Bayonet design
- Lengths of 20, 22, 24, 26cm
- Straight or curved tips

Ultra Fine Micro Scissors

- Extremely thin blades, gold springs
- Ideal for dissection of arachnoid
- Round or flat handle
- Bayonet design
- Lengths of 20, 22, 24, 26cm
- Straight or curved tips

Meningioma Micro Scissors

- Longer blades with a blunt tip (Metzenbaum Style)
- Designed for blunt and sharp dissection
- Round handle
- Bayonet design
- 20cm length
- Straight or curved tips

Meningioma/Tentorium Micro Scissors

- Designed for opening the tentorium
- Short blades with sharp tips, black springs
- Lengths of 20, 22cm
- Bayonet design
- Flat handle
- Straight or curved tips

Lawton Style Micro Scissors

- Round or flat handle
- 18cm Length
- Straight or curved tips
- Non bayonet

Lengths and Tips



65mm working length,
20cm overall



85mm working length,
22cm overall



110mm working
length, 24cm overall



135mm working
length, 26cm overall

Straight Tip



Ultra Fine Straight Tip



Curved Tip



Ultra Fine Curved Tip



Straight Blunt Tip



Curved Blunt Tip



Mizuho America, Inc.
133 Brimbal Avenue
Beverly, MA 01915
www.mizuho.com


Call Toll-free:: 800 699 2547
email:: mizuho@mizuho.com



Meningioma Micro Scissor



Ultra Fine Micro Scissor



Lawton Micro Scissor

Ordering Information

Round Handled Micro Scissor Set

Catalog No.	Description
07-827-01	Micro Scissors Straight 20x65
07-827-02	Micro Scissors Curved 20x65
07-827-03	Micro Scissors Straight 22x85
07-827-04	Micro Scissors Curved 22x85
07-827-05	Micro Scissors Straight 24x110
07-827-06	Micro Scissors Curved 24x110
07-827-07	Micro Scissors Straight 26x135
07-827-08	Micro Scissors Curved 26x135
07-827-00	Micro Scissors Tray

Flat Handled Micro Scissors Set

Catalog No.	Description
07-827-01F	Flat Handled Scissors 20cm Straight
07-827-02F	Flat Handled Scissors 20cm Curved
07-827-03F	Flat Handled Scissors 22cm Straight
07-827-04F	Flat Handled Scissors 22cm Curved
07-827-05F	Flat Handled Scissors 24cm Straight
07-827-06F	Flat Handled Scissors 24cm Curved
07-827-07F	Flat Handled Scissors 26cm Straight
07-827-08F	Flat Handled Scissors 26cm Curved
07-827-00	Micro Scissors Tray

Ultra-Fine Round Handled Micro Scissors Set

Catalog No.	Description
07-827-22	Micro Scissors 20cm Straight Ultrafine Tip
07-827-23	Micro Scissors 20cm Curved Ultrafine Tip
07-827-24	Micro Scissors 22cm Straight Ultrafine Tip
07-827-25	Micro Scissors 22cm Curved Ultrafine Tip
07-827-26	Micro Scissors 24cm Straight Ultrafine Tip
07-827-27	Micro Scissors 24cm Curved Ultrafine Tip
07-827-28	Micro Scissors 26cm Straight Ultrafine Tip
07-827-29	Micro Scissors 26cm Curved Ultrafine Tip
07-827-00	Micro Scissors Tray

Ultra-Fine Flat Handled Micro Scissors Set

Catalog No.	Description
07-827-22F	Micro Scissors 20cm Straight Ultrafine Tip
07-827-23F	Micro Scissors 20cm Curved Ultrafine Tip
07-827-24F	Micro Scissors 22cm Straight Ultrafine Tip
07-827-25F	Micro Scissors 22cm Curved Ultrafine Tip
07-827-26F	Micro Scissors 24cm Straight Ultrafine Tip
07-827-27F	Micro Scissors 24cm Curved Ultrafine Tip
07-827-28F	Micro Scissors 26cm Straight Ultrafine Tip
07-827-29F	Micro Scissors 26cm Curved Ultrafine Tip
07-827-00	Micro Scissors Tray

Additional Scissors Available

Catalog No.	Description
07-134-10	Pistol Grip Micro Scissors, Straight
07-134-11	Pistol Grip Micro Scissors, Curved right
07-134-12	Pistol Grip Micro Scissors Curved left
07-134-13	Pistol Grip Micro Scissors Curved up
07-827-09	Blunt Tip Scissors 20cm Straight
07-827-10	Blunt Tip Meningioma Scissors 20 CM Curved
07-827-15	Micro Scissors 15cm Straight Flat Handle
07-827-18F	Flat Handled 18cm Straight Scissors
07-827-19F	Flat Handled 18cm Curved Scissors
07-827-20	Micro Scissors 15cm Straight Ultrafine Tip
07-827-20F	Micro Scissors 15cm Straight Ultrafine Tip Flat Handle
07-827-21	Micro Scissors 15cm Curved Ultrafine Tip
07-827-21F	Micro Scissors 15cm Curved Ultrafine Tip Flat Handle
09-071-20F	Meningioma Scissors 20cm Straight Flat Handle
09-071-21F	Meningioma Scissors 20cm Curved Flat Handle
09-071-22F	Meningioma Scissors 22cm Straight Flat Handle
09-071-23F	Meningioma Scissors 22cm Curved Flat Handle
07-827-30	Lawton Scissors 18cm Straight
07-827-30F	Lawton Scissors 18cm Straight Flat Handle
07-827-31	Lawton Scissors 18cm Curved
07-827-31F	Lawton Scissors 18cm Curved Flat Handle
07-827-33	Micro Scissor 10cm Straight Shf 60 deg ang Up UF Blade FP
07-827-34	Stevens Tenotomy Scissor Straight 11cm
07-827-35	Stevens Tenotomy Scissor Curved 11cm
07-827-49	Micro Sylvian Fissure Scissor 30 deg 3cm WL Irrigating
07-827-50	Micro Arteriotomy Scissor 60 deg 6cm WL 0.9mm diameter
07-827-51	Micro Arteriotomy Scissor 30 deg 3cm WL 0.9mm diameter
07-827-52	Micro Arteriotomy Scissor 60 deg 3cm WL 0.9mm diameter

Custom Scissors Available!