



MIZUHO[®]


Mizuho[®] Periosteal and Dural Elevators

- Well balanced
- Various angles
- Lightweight handle
- Greater tactile feel

Mizuho® Periosteal and Dural Elevators


Fukushima designed instruments for compressing dura and for separating the periosteum from the bone.

Handles for instruments are lighter and larger for excellent tactile feel.




07-825-01

Langenbeck Elevator, 13mm Width, 20cm Overall Length




07-825-03

Slight Curve Elevator, 6mm Width, 17.8cm Overall Length




07-825-02

Very Slight Curve Elevator, 6mm Width, 17.8cm Overall Length




07-825-04

Full Curve Elevator, 6mm Width, 17.8cm Overall Length




07-826-01

Double Ended Gouge



07-826-03

Reduced Double Ended Elevator and Dissector





07-826-02

Double Ended Elevator and Dissector


Mizuho® Bayonet Probe with Touch Spring

Monopolar cautery specifically indicated for use in transsphenoidal exposure and other skull base procedures.



07-824-01

55mm Working Length



07-824-02

110mm Working Length

Japan
Mizuho Medical Co.,Ltd.
Tel: 81-3-3815-7101

USA
Mizuho America, Inc.
Tel: 1-800-699-2547

Europe
Medical Products International
Tel: 31-(0)-318642834

MIZUHO®